

# NicAlert™

## READING THE DEVELOPED STRIP

**LEVEL 0** 0-10 ng/ml



**LEVEL 1** 10-30 ng/ml



**LEVEL 2** 30-100 ng/ml



**LEVEL 3** 100-200 ng/ml



**LEVEL 4** 200-500 ng/ml



**LEVEL 5** 500-1000 ng/ml



**LEVEL 6** 1000+ ng/ml



**IMPORTANT:** LOWEST PINK BAND READING IS THE FINAL RESULT  
(NOTE: MAY NOT NECESSARILY BE THE DARKEST BAND)